

## ADDING/REMOVING NEW SCREENS

If you want to add a new screen (called component), please follow the steps below:

1. Identify the kind of component you want to add. We have two kinds of components. The first type does not have a provision to add files. It only supports readings. The second type supports both readings and files. Refer to the snapshots below:

Type - I component

|          |                |            |     |              |                                        |                 |                  |
|----------|----------------|------------|-----|--------------|----------------------------------------|-----------------|------------------|
| Metadata | Film thickness | UV-Vis-NIR | I-V | Profilometry | Conductivity (from I-V & profilometry) | <b>Mobility</b> | GIWAXS (NSLS II) |
|----------|----------------|------------|-----|--------------|----------------------------------------|-----------------|------------------|

---

|          |                                                        |         |                                                   |              |                                                           |
|----------|--------------------------------------------------------|---------|---------------------------------------------------|--------------|-----------------------------------------------------------|
| Batch ID | <input type="text" value="Please enter batch number"/> | Reading | <input type="text" value="Please enter reading"/> | Entry Person | <input type="text" value="Please enter your group name"/> |
|----------|--------------------------------------------------------|---------|---------------------------------------------------|--------------|-----------------------------------------------------------|

Type - II component

|          |                |                   |     |              |                                        |
|----------|----------------|-------------------|-----|--------------|----------------------------------------|
| Metadata | Film thickness | <b>UV-Vis-NIR</b> | I-V | Profilometry | Conductivity (from I-V & profilometry) |
|----------|----------------|-------------------|-----|--------------|----------------------------------------|

---

|         |                                                        |              |                                                   |              |                                           |
|---------|--------------------------------------------------------|--------------|---------------------------------------------------|--------------|-------------------------------------------|
| BatchID | <input type="text" value="Please enter batch number"/> | Reading (cm) | <input type="text" value="Please enter reading"/> | Entry Person | <input type="text" value="Please enter"/> |
|---------|--------------------------------------------------------|--------------|---------------------------------------------------|--------------|-------------------------------------------|

---

**Upload**

Drag 'n' drop so

**Files**

**Uploaded Files**

|                                                                                           |             |
|-------------------------------------------------------------------------------------------|-------------|
| <a href="#">0ea402b6-4ba4-46db-80df-90b608df3e11_Spring24_Preanneal_Sample128 (1).csv</a> | <b>Plot</b> |
|-------------------------------------------------------------------------------------------|-------------|

---

|                                                                                          |             |
|------------------------------------------------------------------------------------------|-------------|
| <a href="#">8ca8add7-3477-4762-8471-6f26fc82dd76_Spring24_Postdope_Sample128 (1).csv</a> | <b>Plot</b> |
|------------------------------------------------------------------------------------------|-------------|

---

|                                                                                            |             |
|--------------------------------------------------------------------------------------------|-------------|
| <a href="#">4f5e9626-96b6-4e66-94b2-62f7b529d9ae_Spring24_Postanneal_Sample128 (1).csv</a> | <b>Plot</b> |
|--------------------------------------------------------------------------------------------|-------------|

2. Front End Changes required to add a component. Go to the [file](#)

```
/dope-frontend/src/dialogs/Form.dialog.js
```

Search for `<TabPanel>`. If you want to copy a Type-1 component look for

```
<TabPanel header="Mobility">
```

You need to copy the entire fragment and make changes to all the names as necessary. Let's say this new component is called "X". Your new `<TabPanel>` will look as follows.

```
<TabPanel header="X">
  <div>
    {formData?.X.map((singleX, idx) => (
      <div className="container" key={idx}>
        <div className="field">
          <span className="text">Batch ID</span>
          <textarea className="textbox" value={formData?.X[idx]?.batch}
onChange={(e) => handleXChange(e,idx,'batch')} rows={1} cols={30}/>

          <span className="text">Reading</span>
          <textarea className="textbox" value={formData?.X[idx]?.reading}
onChange={(e) => handleXChange(e,idx,'reading')} rows={1} cols={30}/>

          <span className="text">Entry Person</span>
          <textarea className="textbox" value={formData?.X[idx]?.who}
onChange={(e) => handleXChange(e,idx,'who')} rows={1} cols={30}/>
        </div>
      </div>
    )
  )}

  <div className="button-container">
    <div className="butt">
      <Button label="Add a new X reading" onClick={addX} />
    </div>
    <div className="butt">
      <Button className="butt" label="Save Data"
onClick={()=>saveProduct(allTypes[0])} />
    </div>
  </div>
</TabPanel>
```

3. In the same file, you need to add two more functions at the top of the same file i.e. /dope-frontend/src/dialogs/Form.dialog.js

```
const handleXChange = (e, index, field) => {
  const newX = [...formData.X];
  if(field === 'batch'){
    newX[index].batch = e.target.value;
  }

  if(field === 'reading'){
    newX[index].reading = e.target.value;
  }

  if(field === 'who'){
    newX[index].who = e.target.value;
  }

  setFormData({ ...formData, X: newX });
};

// Add a new mobility field
const addX = () => {
  setFormData({
    ...formData,
    X: [...formData.X, { batch: "", reading: "", who: "" }],
  });
};
```

If you want to remove a component, say Y, you just need to comment out the associated tab panel for Y e.g.

```
<TabPanel header="Y">
.....
.....
</TabPanel>
```